

Package ‘rawReader’

April 11, 2016

Type Package

Title Read Thermo .RAW files in R

Version 0.1

Date 2016-02-24

Author Tom Wilson <tpw2@aber.ac.uk>

Maintainer Tom Wilson <tpw2@aber.ac.uk>

Description A series of C++ functions are used to access the Thermo MSFilerReader.dll. Each C++ function is designed for single file and single scan input; corresponding R functions and system calls are used as multi scan and multi file wrappers to the C++ functions, via the pre-compiled binary executables.

License GPL (>= 3)

URL <https://github.com/wilsontom/rawReader>

BugReports <https://github.com/wilsontom/rawReader/issues>

LazyData TRUE

RoxygenNote 5.0.1

R topics documented:

compileAll	2
getBValues	2
getCValues	3
getITime	4
getMzInt	5
getMzIntNoise	5
getResComp	6
getRFComp	7
getSpaceComp	7
testPackage	8
Index	9

compileAll	<i>Compile .cpp files to executables (.exe)</i>
------------	---

Description

Wrapper function which compiles all .cpp files into executables (.exe). The Microsoft Visual Studio vcvars32.bat file is used to set include paths and environment variables, so that this function can be simply executed in the Windows cmd line from an R session, as opposed to manually compiling all .cpp files through the Visual Studio C++ Developers Command Prompt. See <https://github.com/wilsonsom/rawReader/blob/master/README.md> for further details on installation.

Usage

```
compileAll(filepath)
```

Arguments

filepath the filepath of the rawReader source directory.

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:
options(VSPATH = "<path_to_visual_studio_vcvars32.bat>")
source("rawReader/R/compileAll.R")
compileAll("C:/rawReader")

## End(Not run)
```

getBValues	<i>Get B Value</i>
------------	--------------------

Description

Get the FT-MS conversion parameter B for each scan in a given range

Usage

```
getBValues(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a numeric vector of the B value for each scan

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getBValues(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

getCValues	<i>Get C Value</i>
------------	--------------------

Description

Get the FT-MS conversion parameter C for each scan in a given range

Usage

```
getCValues(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a numeric vector of the C value for each scan

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getCValues(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

getITtime*Get Ion Injection Time*

Description

Get the ion injection time (IT) for each scan in a given range

Usage

```
getITtime(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a numeric vector of the ion injection time for each scan

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getITtime(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

getMzInt	<i>Get Mass and Intensity</i>
----------	-------------------------------

Description

Get the profile data (m/z and intensity) across a given scan range

Usage

```
getMzInt(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a list of matrices for m/z and intensity.

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getMzInt(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

getMzIntNoise	<i>Get Mass, Intensity and Noise</i>
---------------	--------------------------------------

Description

Get the profile data (m/z, intensity and noise) across a given scan range

Usage

```
getMzIntNoise(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a list of matrices for m/z, intensity and noise

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getMzIntNoise(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

getResComp

Get Resolution Compensation

Description

Get the Resolution mass compensation (ppm) for each scan in a given range

Usage

```
getResComp(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a numeric vector of the resolution compensation for each scan

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getResComp(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

`getRFComp`*Get RF Compensation*

Description

Get the RF mass compensation (ppm) for each scan in a given range

Usage

```
getRFComp(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a numeric vector of the RF compensation for each scan

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getRFComp(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

`getSpaceComp`*Get Space Charge Compensation*

Description

Get the space charge mass compensation (ppm) for each scan in a given range

Usage

```
getSpaceComp(filename, scans = c())
```

Arguments

filename	a .raw file
scans	a numeric vector of scan numbers to extract

Value

a numeric vector of the space charge compensation for each scan

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
getSpaceComp(QC.raw, scans = c(2:24))  
  
## End(Not run)
```

testPackage

Package Testing

Description

Run some basic tests to ensure there have been no problems during compilation

Usage

```
testPackage(scans = c(1:5))
```

Arguments

scans a numeric vector of scan numbers to extract

Author(s)

Tom Wilson <tpw2@aber.ac.uk>

Examples

```
## Not run:  
testPackage(scans = c(1:5))  
  
## End(Not run)
```


Index

`compileAll`, [2](#)

`getBValues`, [2](#)

`getCValues`, [3](#)

`getITtime`, [4](#)

`getMzInt`, [5](#)

`getMzIntNoise`, [5](#)

`getResComp`, [6](#)

`getRFComp`, [7](#)

`getSpaceComp`, [7](#)

`testPackage`, [8](#)